



YOUR SAFETY IS OUR PRIORITY

Handel Street, your number one local PPE supplier

SAHPRA
SOUTH AFRICAN
HEALTH PRODUCTS
REGULATORY AUTHORITY

SAMED
member and proudly committed
to the Medical Device Code



Handel Street
Medical

MDM
MEDICAL DEVICE
MANUFACTURERS
ASSOCIATION
OF SOUTH AFRICA

**SAHPRA
COMPLIANT**

ABOUT US

At Handel Street we focus on providing high quality protective gear to customers across South Africa. We offer a diverse range of products, suited to some of the most extreme work environments.

We strive to fulfil any order, no matter how big or how small, with nothing but the very best in top quality, regulated PPE. If you're looking for quality and good value for money, then we have got you covered. Give us a call to see how we can help you.

OUR SERVICES



Medical PPE Equipment

Handel Street supplies quality personal protective equipment to the South African Medical industry.



Industrial PPE Equipment

Handel Street supplies top rated personal protective equipment to the South African Industrial industry.



**PROUDLY PRODUCED
IN SOUTH AFRICA**



GIVE US A CALL ON:
+27 (0)31 332 6029



37 Cato Street, Durban, Kwa-Zulu Natal

Email: hsdurban@handelstreet.co.za




www.shophandelstreet.co.za



Level 1A Mask 4 PLY[®] PATIENT & HCW SURGICAL MASK

Handel Street, your number one local PPE supplier

 **PROUDLY PRODUCED IN SOUTH AFRICA**



PRODUCT INFO

Model : FACE FIT LVL1 GC MASK

Part Codes Applicable:

HSBFESM3 NAPPI CODE APPLICABLE : 1068891001 (SMALL)

HSBFESM4 NAPPI CODE APPLICABLE : 1068912001 (MEDIUM)

HSBFESM5 NAPPI CODE APPLICABLE : 1068912001 (LARGE)

Surgical Disposable Face Mask

95% BFE protection 95% PFE

GMDN Code 35177

Efficient against aerosols in the ranges of 0.1

PRODUCT FEATURES

Face seal : Yes

Nose foam : Yes

FDA Cleared : Yes

EN14683 Compliant : Yes

SANS1866-1 : 2018 Compliant : Yes

Features Advanced Electrostatic Media

Natural Rubber Latex Components: No

Number per Box: 25 Individually Wrapped

Packaging Bulk : 450 Units Per Carton 18 Boxes Per Carton

Product Color : White

Product Type : Healthcare Surgical Mask

FLAME SPREAD CLASS : 16 CFR

DUNS Number 56-887-6648

 Call us on: +27 (0)31 332 6029

 37 Cato Street, Durban, Kwa-Zulu Natal

 Email: hsdurban@handelstreet.co.za

 www.shophandelstreet.co.za



YOUR SAFETY IS OUR PRIORITY

Handel Street, your number one local PPE supplier

 **PROUDLY PRODUCED IN SOUTH AFRICA**

SAHPRA
SOUTH AFRICAN
HEALTH PRODUCTS
REGULATORY AUTHORITY

SAMED
member and proudly committed
to the Medical Device Code


Handel Street
Medical

MDM MEDICAL DEVICE
MANUFACTURERS
ASSOCIATION
OF SOUTH AFRICA

FITTING INSTRUCTIONS

Acute Respiratory Sicknesses can be transmitted by respiratory droplets over a short distance or through direct contact with a patient's secretions. Wearing a mask offers protection against various airborne bacteria, germs and viruses. If you have a respiratory tract infection, it also helps prevent the spread of the illness. Surgical mask, if properly worn, are effective in preventing the spread of droplet infections. Wearing a mask is just one way to help prevent respiratory tract infections. Most important is to observe good personal hygiene. Wash hands before putting on a mask, before and after taking one off. When wearing a surgical mask, ensure that: The mask fits snugly over the face. The front side of the mask faces outwards. The strings or elastic bands are positioned properly to keep the mask firmly in place. The mask covers the nose, mouth and chin. Try not to touch the mask once it is secured on your face as frequent handling may reduce its protection. If you must do so, wash your hands before and after touching the mask. When taking off the mask, avoid touching the outside of the mask as this part may be covered with germs. After taking off the mask, fold the mask outwards (i.e. the outside of the mask facing inwards), then put the mask into a plastic or paper bag before putting it into a rubbish bin with a lid. A surgical mask should be discarded after use and under no circumstances should it be used for longer than a day. Replace the mask immediately if it is damaged or soiled.



GIVE US A CALL ON:
+27 (0)31 332 6029

 37 Cato Street, Durban, Kwa-Zulu Natal

 Email: hsdurban@handelstreet.co.za

 www.shophandelstreet.co.za